

## ZEISS Teleloupe Spectacles

The magnifying system for professional applications



Teleloupe Spectacles GTX

### Teleloupe Spectacles GTX

The Teleloupe Spectacles GTX feature some outstanding benefits: They are fitted almost like a normal pair of glasses by your eyecare professional. Both short-sightedness (myopia) and far-sightedness (hyperopia) can be corrected. Objects in the close-up range are seen stereoscopically with 2x or 2.5x magnification.



Teleloupe Spectacles GTX with Saphiro<sup>2</sup>

The Teleloupe Spectacles GTX are available with several fixed working distances between 300 and 500 mm (2x magnification) and 350 mm to 450 mm (2.5x magnification). The Teleloupe Spectacles GTX can be equipped with the Saphiro<sup>2</sup> LED illumination. This allows an ideal, very easy combination of magnification and illumination.



Teleloupe Spectacles G 2.5 TTL with sportsframe

## Teleloupe Spectacles G 2.5 TTL

### Working with precision down to even the finest detail

Outstanding image quality over a large field of view is the typical feature of all magnifying visual devices from ZEISS.

Teleloupe Spectacles G 2.5 TTL offers a field of view which is up to 30% larger than the GTX system.

### Technical data:

- Magnification 2.5x
- Working distances: 350 / 400 / 450 mm
- Fields of view: 75 / 85 / 95 mm
- Individually matched to specific requirements
- Relaxed work without strain
- Crafted to individual centring data

### Availability

- **System carrier:** STMS Titanium frame (with straight temples and headband)  
Sportsframe (plastic frame with straight temples and headband; curved lenses provide a splash guard)
- **Individual prescription:** STMS can be glazed with prescription lenses  
Sportsframe comes with a clip-on for prescription lenses



Teleloupe Spectacles G 2.5 TTL with Saphiro<sup>2</sup>

The G 2.5 TTL Teleloupe Spectacles can be combined with Saphiro<sup>2</sup> LED illumination. This results in an ideal combination of magnification and illumination.

ZEISS is one of the world's leading manufacturers of spectacle lenses, and is committed to delivering maximum precision and comfort. ZEISS designs and produces lenses, instruments and measurement systems, as well as retail concepts and technology services that continue to raise the bar in vision care.

### Contact us to get started

ZEISS Customer Service

+27 11 538 4200

Email





---

**Explore**

- [Understanding Vision](#)
- [Health + Prevention](#)
- [Lifestyle + Fashion](#)
- [Driving + Mobility](#)
- [Sports + Leisure](#)
- [Work Life](#)

**Help me choose**

- [Distance + Reading Glasses](#)
- [Varifocal Lenses](#)
- [Sunglasses](#)
- [Working Glasses](#)
- [Sports Glasses](#)
- [Glasses for Children](#)
- [Lens Coatings](#)
- [Lens Cleaning](#)
- [At the Optician](#)

**Services**

- [My Vision Profile](#)
- [Online Vision Check](#)

**For Eye Care Professionals**

- [Instruments + Technologies](#)
- [ZEISS Spectacle Lenses](#)
- [ZEISS Cleaning Solutions](#)